



# MONARK 927X



## MAIN FEATURES

- User friendly exercise bike for rehabilitation
- Easy entry frame
- Extra durable design

## THE BASICS

- Good cycle feeling
- 5 kHz HR systems
- User friendly display
- Display with 12 training protocols  
HR control and user programs
- HR in handlebar

Max user weight 250 kg

Medical Device class IIa

# Specifications

## MONARK 927X

<b>Certifications</b>	EU directive 93/42/EU class IIa
<b>Frame</b>	
<b>Type</b>	Medical Frame
<b>Handlebar</b>	Monark handlebar with pulse grips
<b>Stem</b>	n/a
<b>Seatpost</b>	Monark Standard
<b>Saddle</b>	Saddle, black
<b>BB</b>	68/122 mm
<b>Crank</b>	Steel, 170 mm
<b>Pedals</b>	9/16", pedals with foot straps
<b>RESISTANCE</b>	
<b>System</b>	Spring
<b>Flywheel</b>	20 kg
<b>DISPLAY</b>	
<b>RPM</b>	Yes
<b>HR</b>	Handlebar with pulse grips
<b>Time</b>	Yes
<b>Speed</b>	Yes
<b>Distance</b>	Yes
<b>Watt</b>	Yes
<b>Connectivity</b>	
<b>RS232</b>	n/a
<b>Usb</b>	n/a
<b>WEIGHT &amp; MEASUREMENTS</b>	
<b>Max user weight</b>	250 kg
<b>Lenght (mm)</b>	1240 mm
<b>Width (mm)</b>	500 mm at handlebar, 640 mm at support legs
<b>Height (mm)</b>	1260 mm at handlebar 780-1170 mm at saddle
<b>Weight (kg)</b>	56 kg
<b>POWER REQUIREMENTS</b>	220V (9V, 500 mA)



### HEART RATE

Equipped to measure heartrate via chestbelt or sensors in the handles.

### MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 153 601 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden