




Belastningsanordning  
höger sida



Belastningsanordning  
vänster sida



Pos.	Qty.	Art. No.	Description	Pos.	Qty.	Art. No.	Description
1	2	9125-86	Spring	9	1	9384-27	Pulley
2	1	9311-73	Brake belt LT, complete	10	1	9326-164	Magnet
3	1	9300-3	Flywheel, complete	11	1	9300-475	Crank set, 52T, 172,5 mm, Q 146 mm
4	1	9300-24	-Wheel suspension complete set	12	1	9300-480	Cartridge bottom bracket 68/110 mm
5	1	9100-20	Spring for press arm	13	1	9310-93	Inner chain guard w sensor and nuts
6	1	9311-59	Press arm	14	1	9311-161	-Sensor with cable
7	1	9328-43	Tension lever	15	1	9310-118	Chain 9300, 120 L, with chain lock
8	1	9311-67	Potentiometer				






Fig: Chest belt with pulse transmitter (accessory)



Fig: Power adaptor SE (accessory)

### LT2 white

Pos.	Qty.	Art. No.	Description	Pos.	Qty.	Art. No.	Description
1	1	9334-110	Saddle, blue		1	9311-3	Handlebar with stem and clamp, complete
	1	9336-2	Saddle adaptor 22 mm	12	1	9311-33	-Handlebar stem
2	1	9311-23	Sledge, complete	13	1	9311-31	-Handlebar clamp
3	1	9311-24	-Locking handle M10x32, black	14	1	9311-24	-Locking handle M10x32, black
4	1	9311-21	Saddle post	15	1	C2305771-42	-Reparto corse handlebar
5	2	8321-75	Pedals SPD with clips and strap		1	C2600206	-Handlebar tape, black, CORK
6	1	9301-3	Support tube rear, white, complete	16	1	9311-4	Workload control unit, complete
	2	9328-51	-Plastic cap, blue		1	9311-44	-Wire adjustment
	2	9328-26	-Rubber foot with screw M8		1	9311-45	-Rubber ball 32/M8, black
7	1	9301-4	Support tube front, white, complete	17	1	9311-164	Display LT2
	2	9328-51	-Plastic cap, blue	18	2	9000-104	Pole screw M5x12, black
	2	9328-26	-Rubber foot with screw M8	19	1	9326-801	End cap with hole, black
	2	9328-37	-Transport wheel, complete	20	1	9311-610	Cover, white
8	1	9310-595	Side cover rear, right, white	21	1	9310-605	Side cover rear, left
9	2	9000-103	Pole screw M5x12, white	22	1	9000-103	-Pole screw M5x12, white
10	2	9310-27	Piston locking, white, complete		1	9384-650	Power adaptor SE (accessory)
	2	9100-289	-Locking handle blue, M10x50		1	9384-62	Power adaptor other countries (accessory)
11	2	9328-315	Bushing with hole, blue		1	9311-75	Chest belt (accessory)
					1	9000-211	Calibration weight 4 kg (accessory)
					1	9000-212	Calibration weight 2 kg (accessory)

### LT2 black

Same as above except for a few positions that differ

Pos.	Qty.	Art. No.	Description	Pos.	Qty.	Art. No.	Description
1	1	9336-115	Saddle, black, complete	8	1	9310-597	Side cover rear, right, black
6	1	9301-167	Support tube rear, black, complete	9	2	9000-104	Pole screw M5x12, black
	2	9309-51	-Plastic cap, black	10	2	9310-275	Piston locking, black, complete
	2	9328-26	-Rubber foot with screw M8		2	9100-287	-Locking handle M10x50, black
7	1	9301-157	Support tube front, black, complete	11	2	9309-1315	Bushing with hole, black
	2	9309-51	-Plastic cap, black	20	1	9311-167	Cover, black
	2	9328-26	-Rubber foot with screw M8	21	1	9310-607	Side cover rear, left, black
	2	9328-37	-Transport wheel, complete	22	1	9000-104	-Pole screw M5x12, black